

# Specification Sheet

Product Name	<b>CRYO-DMSO-F</b> , liquid, sterile-filtered, 100 mL
Product Number	10005-01
Product Brand	Revive Organtech, Inc
Storage Temp	2-8°C
<b>TEST</b>	<b>SPECIFICATION</b>
Appearance (Turbidity)	Clear liquid
Appearance (Form)	Solution
pH	7.80~8.20
Sterility	Sterile Filtered (0.22uM)
Endotoxin Level	N/A
<b>Instructions for use</b>	<p>CRYO-DMSO-F is a serum-free, protein free, animal origin-free, <b>DMSO Free</b> and fully defined cryopreservation medium. Designed to prepare and preserve cells in ultralow temperature environments (-80 to -196°C), CRYO-DMSO-F provides a safe, protective environment for cells during the freezing, storage, and thawing process. CRYO-DMSO-F contained proprietary components which are directly reducing the level of freezing induced apoptosis and necrosis and improving post-thaw cell viability and function. CRYO-DMSO-F is recommended for the preservation of Primary cells, iPSC, MS, MSC, Oocyte, Embryo, CAR-T, NK, and Macrophage cells. CRYO-DMSO-F is cGMP-mimic condition manufactured with high quality grade components.</p> <p><b>Preparation instructions</b></p> <p>The CRYO-DMSO-F solution is ready-to-use and complete with no additives required. Wipe down the outside of container with 70% alcohol before opening as the contents are sterile. If the seal has been broken, do not use it.</p> <p><b>Storage/Stability</b></p> <p>Store the CRYO-DMSO-F solution at 2-8°C and protected from light until ready to use it.</p> <p><b>Freezing Procedures</b></p> <ol style="list-style-type: none"> <li>1. Suspended cell to be cryopreserved using mechanically or enzymatically dissociation.</li> </ol>

2. Centrifuge the cells to obtained pellet.
3. Remove supernatant (remove the culture media as possible to reduce dilution of the CRYO-DMSO-F solution).
4. Add Ambient temperature CRYO-DMSO-F solution to a cell concentration range of  $1-10 \times 10^7$  cells / 1 ml of CRYO-DMSO-F for standard cell culture protocol.
5. After mixed with CRYO-DMSO-F solution with cells, incubated for 10 minutes at  $1-4^{\circ}\text{C}$  to penetrate cryoprotectants inside of **cells (in case of small tissue or organoid, 20 minutes incubation at  $1-4^{\circ}\text{C}$ )**
6. Nucleation-lower sample temperature  $-80^{\circ}\text{C}$ ; After cells are mixed with solution, put cryovial into controlled rate freezer ( $-1^{\circ}\text{C}/\text{minute}$ ) and then freeze to  $-80^{\circ}\text{C}$  (slow freezing method), or put cryovial into Bicell, Mr. Frosty Freezing container, or similar kinds of slow freezing container and put the cryovial included in such slow freezing container into  $-80^{\circ}\text{C}$  freezer.

#### **Storage**

1. After finished the nucleation of cells at  $-80^{\circ}\text{C}$  freezer or controlled rate freezer, store the freeze cell at liquid nitrogen tank (below  $-130^{\circ}\text{C}$ ).

#### **Thawing Procedures**

1. Thaw samples quickly in a  $37^{\circ}\text{C}$  water bath. Samples should be thawed with gentle swirling of the sample until all visible ice has melted (*Do not allow sample to warm above chilled temperatures ( $0-10^{\circ}\text{C}$ ). Cryovials should be cool to the touch when removed from the water bath.*)
2. Dilute cell/CRYO-DMSO-F mixture immediately with appropriate culture medium. Add 20 ml of culture medium to 50 ml conical tube and gently mix with thawed cell/CRYO-DMSO-F mixture to culture medium (The dilution culture medium should be  $20-37^{\circ}\text{C}$ ). **A dilution ratio of 1:20 V: V (1 mL of thaw cell + 20 ml of culture media)** or greater is recommended. After diluting cells with warmed culture medium, gently inverted for 5-10 times and follow by centrifuge (recommend at 400G for 5 minutes).
3. After centrifuge, completely suction of supernatant and add new warmed culture media.
4. Plate the cells appropriately and cultured the cells or use immediately.